

# NexScope

NEPHROSCOPE AND MINI NEPHROSCOPE




**NEW**

COMPETENCE IN SURGICAL LASER SINCE 1989  
LASER TECHNOLOGY - MADE IN GERMANY  
[WWW.LISALASER.DE](http://WWW.LISALASER.DE)



## KEY FEATURES

- ★ Outstanding performance: contrast enhanced image bundle, best vision
- ★ Safe and reliable: robust, tightly sealed, laser welded, easy to use and to clean
- ★ Direction of view ideally aligned at 8°
- ★ Strong irrigation flow enables high stone-free rates
- ★ Long lifetime for lowest cost of ownership
- ★ Lightweight design
- ★ Detachable instrument bridge for optimized cleaning and inspection
- ★ Compatible with most popular laser and lithotripter systems
- ★ Autoclavable: 18 min. @ 134°
- ★ Made in Germany



## Technical Data **NexScope**

### Nephroscope and Mini Nephroscope including detachable instrument bridge

Nephrosopes	Working channel	Field of View	Part number
Mini NexScope 13.5 F 220	6 F	84°	101 503 435
Mini NexScope 13.5 F 320	6 F	84°	101 503 436
NexScope 19 F 220	10 F	84°	101 503 437

Accessories	Description	
NexScope Instrument Bridge	for mini nephroscope and nephroscope	101 503 438
NexScope Sheath 18 F	with obturator for mini nephroscope	101 503 439
NexScope Sheath 24 F	with obturator for nephroscope	101 503 440
NexScope Cleaning Brush For Working Channel	blue, for nephroscope	101 503 486
NexScope Cleaning Brush For Working Channel	white, for mini nephroscope	101 503 487
NexScope Cleaning Brush For Side Port	white, for instrument bridge side port	101 503 488
NexScope Cleaning Brush For Main Port	blue, for instrument bridge main port	101 503 489
Duck Bill Valves	disposable, package á 10	101 503 498
Sealing Caps	disposable, package á 10	101 503 499

Sheath and Obturator



Instrument Bridge



Cleaning Brush



CE acc. Council Directive 93/42/EEC

U.S. federal law restricts these devices to sale by or on the order of a physician

Specifications are subject to change without notice.  
Made in Germany 2015 - 01  
Brochure Part No. 036 013 013



LISA laser products OHG  
Max - Planck - Str. 1  
D - 37191 Katlenburg-Lindau  
Deutschland  
fon: +49 (0)5556 - 9938 - 0  
fax: +49 (0)5556 - 9938 - 10  
info@lislaser.de www.lislaser.de



**LISA LASER PRODUCTS**  
COMPETENCE IN SURGICAL LASER

